PHARMACY STORAGE

BREAK OUT/RECEIVING

Shipments from drug companies are unpacked in the Break Out/Receiving area. Modular Casework cabinets, shelving units and workstations allow staff as well as incoming inventory to become and remain organized.

STAFF LOUNGE

The Staff Lounge is a great place to take a break. It's also a great place to store and secure your personal stuff in lockable laminate or steel lockers.

NARCOTICS

The Narcotics area is where all controlled substances are stored. Lockable cabinets and drawers are critical to ensuring these substances stay out of the wrong hands.

CHECKING

Pharmacists check medication orders before they are sent out in the Checking area. Standing height workstations provide ergonomic and organizational benefits and help ensure that the right medication gets to the right patient.

IV STORAGE

IV Fluids are kept in the IV Storage area. Wire baskets on high-density mobile or FrameWRX shelving help keep those IV's organized and accessible.

CHEMO/COMPOUNDING

The mixing process for topicals, ointments, creams and other liquids takes place in the Chemo/Compounding area. This process can be tough on furniture finishes. Modular Casework with phenolic resin tops stand up to these materials and harsh cleaning fluids.







BULK STORAGE

General and overflow storage happens in the Bulk Storage area. Depending on accessandspaceconstraints, high-density mobile shelving or static shelving will do the trick for storing excess boxes and packaged medication.

PHARMACY WORK STATIONS

From the Pharmacy work area, it is important for the pharmacist to be able to see the entire pharmacy. Seated modular casework workstations and demountable walls allow a clean line of site throughout the space.

UNIT DOSE/ACTIVE STORE

There are several ways to approach storage in the Unit Dose/Active Store/Filling area. Vertical Carousels provide securestorage with quick access. Powered high-density mobile storage with storefront shelving is a great use of space while keeping the most frequently used items close at hand. Modular Casework workstations with half-height medication drawers provide an efficient place to get the job done.

ANTE ROOM/STERILE AREA

Inside Ante Rooms and Sterile areas, contamination must be kept at bay. Count on FrameWRX shelving and wire baskets to keep IV's easy to stock, visible, clean and organized.

PACK AREA

Bulk medications are packed into unit doses in the Pack Area, a job that is often performed by a robot. To assist in this process and to store packaging materials, FrameWRX shelving with EZ Rails and pegs are a good choice.

SOLVING STUFF PROBLEMS FOR PHARMACIES